

# SHIELD<sup>me</sup>™



## DISPOSABLE BOUFFANT CAP

SDM-HCAP



**SHIELD<sup>me</sup>™ NON-STERILE DISPOSABLE BOUFFANT SCRUB CAPS** is also called disposable hair cover, disposable head cover, disposable nurse cap, disposable medical caps, which is the most popular disposable non-woven caps

### SPECIFICATIONS:

- Fabrics** : PP Spun Bond max upto 25gsm  
**Size** : 18 inches, 21 inches, 24 inches with 2 LP Elastic  
**Color** : White/blue/green/yellow/red/black  
**Packing** : 2000 pieces/Carton 54x26x23cm

### COMFORT FITTINGS

- Heat sealing, avoids escape of hair, dandruff and microorganisms during process
- Water-resistant
- Economic



شركة البيع للمعدات واللوازم ش.م.م.  
**Al Nab'a Equipment & Supplies LLC**

Phone: 24229483, Mobile: 99889341

Email: enquiry@alnabaequipment.com, www.alnabaequipment.com

\*We reserve the right to change specifications of our product without prior notice.