

SHIELD^{me}™

CE

DISPOSABLE BOUFFANT HOOD



SHIELD^{me}™ NON-STERILE DISPOSABLE BOUFFANT HOOD The single use SHIELDme Non-Sterile Disposable Hood has design to ensure a perfect fit. Made of 100% polypropylene breathable SMS fabric, which gives the wearer cloth-like comfort. The hood features an extra-long yoke for maximum coverage when worn in conjunction with a coverall with collar, and features an elasticated face-opening with reinforced edges to avoid contamination entering the controlled environment.



SPECIFICATIONS:

- Weight 25-50 g/m²
- Single elasticated front
- V neck and bib
- Light and breathable
- White colour
- Uni-Size
- Comfort fittings

COMFORT FITTINGS

- Resistant to blood and body Fluid Penetration
- Resistant to viral penetration
- Resistant to Blood Borne Pathogens
- Resistant to Chemical, oil, bleach and 50% Sodium Hydroxide



شركة البيع للمعدات واللوازم ش.م.م.
Al Nab'a Equipment & Supplies LLC

Phone: 24229483, Mobile: 99889341

Email: enquiry@alnabaequipment.com, www.alnabaequipment.com

**We reserve the right to change specifications of our product without prior notice.*