

SHIELDme™



NON-COTACT IR THERMOMETER

SME-IRTM01

CE

SHIELDme™ NON-COTACT INFRARED THERMOMETER is a medical device for body / forehead temperature checking, super fast one-second reading with high measurement accuracy.

FEATURES:

With its innovative no-touch technology, this infrared thermometer is suitable for family members of all ages. It has smart built-in fever alarm system.

SPECIFICATION PARAMETER:

Product Name	Non-Contact IR Thermometer
Product Model	SHIELDme
Body Temperature Range	34.0°C ~ 42.9 °C / 93 - 109 °F
Object Temperature Range	0°C - 100°C
Display Resolution	±0.2°C (0.2°F)
Measurement Distance	1 ~ 5 cm
Measurement Time	0.5 s
Automatic Power-off Function	15 s
Relative Humidity	≤85%
Storage Condition	-20°C ~ 50°C
Power Supply	DC 3.0V (2*AAA batteries)
Item Size	116 * 43 * 152mm (L * W * H)



شركة البيع للمعدات واللوازم ش.م.م.
Al Nab'a Equipment & Supplies LLC

Phone: 24229483, Mobile: 99889341

Email: enquiry@alnabaequipment.com, www.alnabaequipment.com

*We reserve the right to change specifications of our product without prior notice.